

Porous Gelatin Methacryloyl

Porous GelMA

Product component

| Item | Character | Package Size | Notes |
|-----------------------|--------------|----------------|--------------|
| Porous GelMA hydrogel | White spongy | 0.6g/bottle x2 | Keep in dark |

This instruction applies to EFL-GM-PR-001/002

Storage

Dry kit: -20℃, 18 months; 4℃, 3 months. (Due to its special properties, storage at -20℃ is recommended). **Sterile solution:** 4℃ (in dark), 7 days; -20℃ (in dark), 6 months. Please note that repeated freezing and thawing of the solution will affect the performance of the product, so it is best to prepare it when using it.

Period of validity

The date of manufacture is shown in the package.

Required materials

- EFL-GM-PR series porous GelMA hydrogel products^{EFL}
- EFL-LS-1601 series 405nm curing light source equipment^{EFL}
- PBS (1X)
- Constant temperature magnetic stirring water baths
- 0.22μm Sterile needle filters
- 10-50mL Sterile centrifuge tubes
- 10mL Syringes
- 1-5mL Pipettes & tips



企业微信公众号
扫描右侧二维码
获取更多信息

Operation steps (3D cell culture)

| Steps | Title | Materials | Processes |
|-------|--------------------|---|--|
| 1 | Prepare solution | <ul style="list-style-type: none"> ➤ EFL-GM-PR series porous GelMA hydrogel products ➤ PBS (1X) ➤ Pipette guns ➤ constant temperature magnetic stirring water baths | <ol style="list-style-type: none"> 1) Add appropriate amount of PBS to the porous GelMA bottle; Recommended concentration for porous GelMA is 6-8% (w/v). (bottle contains magnetic rotor and 0.6g porous GelMA product). 2) Preparation of porous GelMA hydrogel precursor solution by magnetic stirring at 37°C in a water bath keeping in the dark for 1h (important step). |
| 2 | Sterilise solution | <ul style="list-style-type: none"> ➤ 0.22µm Sterile needle filters ➤ Constant temperature water baths ➤ Sterile centrifuge tubes | <ol style="list-style-type: none"> 1) Sterilise the above solution immediately with a sterile 0.22µm needle filter (to prevent gelation when the temperature drops) and store at 37°C in the dark. Note: If it is not possible to use it all at once, store in the refrigerator for a short period of time (< 7 days). Before the next use, redissolve at 37°C and vortex for 20-30 seconds to homogenise the material. |
| 3 | Mix cells | | <ol style="list-style-type: none"> 1) Collect cells. 2) Resuspend cells in the sterile precursor solution (multiple blowing or shaking). |



企业微信公众号
扫描右侧二维码
获取更多信息

| 4 | Cure GelMA | <ul style="list-style-type: none"> ➤ Pipettes & tips ➤ EFL-LS-1601 series 405nm curing light source equipment | <ol style="list-style-type: none"> 1) Add precursor solution to the well plate; 96-well plate: 50-100μL/ well, 48-well plate: 100-300μL/ well, 24-well plate: 300-500μL/ well. 2) Let stand at room temperature for 2 min. 3) Irradiation with EFL-LS-1601 series 405nm light source to cure hydrogels, Irradiation times are shown in the table below: <table border="1"> <thead> <tr> <th>Models</th><th>6%</th><th>7%</th><th>8%</th></tr> </thead> <tbody> <tr> <td>001</td><td rowspan="2">16~18s</td><td rowspan="2">13~16s</td><td rowspan="2">8~10s</td></tr> <tr> <td>002</td></tr> </tbody> </table> | Models | 6% | 7% | 8% | 001 | 16~18s | 13~16s | 8~10s | 002 |
|--------|---------------|---|---|--------|----|----|----|-----|--------|--------|-------|-----|
| Models | 6% | 7% | 8% | | | | | | | | | |
| 001 | 16~18s | 13~16s | 8~10s | | | | | | | | | |
| 002 | | | | | | | | | | | | |
| 5 | Wash samples | <ul style="list-style-type: none"> ➤ Pipettes & tips | <ol style="list-style-type: none"> 1) Add medium and incubate at 37°C for 5 minutes. 2) Remove the medium. | | | | | | | | | |
| 6 | Culture cells | | <ol style="list-style-type: none"> 1) Change medium, observe and photograph according to experimental design. | | | | | | | | | |

Operation steps (2D cell culture)

The main steps of 2D cell culture are the same as those of 3D culture, except that step 3 mixing cells. After washing the sample in step 5, seed cells on the surface.



企业微信公众号
扫描右侧二维码
获取更多信息